

ts/isotope_products/isotrak_calibration_sources/reference_sources/planchet_disk_sources/alpha_sources_af_type/a2_capsule_type/)

Home (home/) / Products (home/products/) / Isotrak Calibration Sources (home/products/isotope_products/isotrak_calibration_sources/) /

Reference Sources (home/products/isotope_products/isotrak_calibration_sources/reference_sources/) /

Planchet & Disk Sources (home/products/isotope_products/isotrak_calibration_sources/reference_sources/planchet_disk_sources/) /

Alpha sources AF-Type (home/products/isotope_products/isotrak_calibration_sources/reference_sources/planchet_disk_sources/alpha_sources_af_type/) /

A2 capsule type (home/products/isotope_products/isotrak_calibration_sources/reference_sources/planchet_disk_sources/alpha_sources_af_type/a2_capsule_type/)


A2 capsule type

Alpha Particle Standards—Type A-2

The type A-2 source is permanently fixed in an aluminum holder 0.5" diameter x 0.250" high (12.7 mm x 6.35 mm). The active diameter is 0.197" (5.0 mm).



(fileadmin/_processed_/csm_catalogue_reference_and_calibration_sources_Seite_37_Bild_0001_0397b926c1.jpg)



(fileadmin/_processed_/csm_alpha_A2_cc616fd905f.jpg)

All alpha standards are offered as spectral grade sources up to the

activity and active diameters listed unless otherwise noted. All electroplated alpha standards are manufactured to a tolerance of +/-30% of the nominal activity.

Nuclide	Product code	Activity range	Remark
Am-241	AF-241-A2	185 Bq - 37 kBq	
Cf-252	AF-252-A2	185 Bq - 37 kBq	

Cm-244	AF-244-A2	185 Bq - 37 kBq	
Gd-148	AF-148-A2	185 Bq - 3,7 kBq	
Np-237	AF-237-A2	185 Bq	maximum activity, Not supplied as spectroscopy grade sources.
Po-210	AF-210-A2	185 Bq - 37 kBq	
Pu-238	AF-238P-A2	185 Bq - 37 kBq	
Pu-239	AF-239-A2	185 Bq - 370 Bq	
Ra-226	AF-226-A2	185 Bq - 3,7 kBq	Not necessarily in equilibrium with daughters at time of shipment.Ra-226 will leak Rn-222 gas.
Th-228	AF-228-A2	185 Bq - 18,5 kBq	Not necessarily in equilibrium with daughters at time of shipment.Th-228 will leak Rn-220 gas.
Th-230	AF-230-A2	185 Bq - 3,7 kBq	
Th-232	AF-232-A2		Not available
U-235	AF-235-A2	15 Bq	maximum activity, Not supplied as spectroscopy grade sources.Not necessarily in equilibrium with daughters at time of shipment.
U-238 (Nat)	AF-238(Nat)-A2	1,85 Bq	maximum activity, Not supplied as spectroscopy grade sources.Not necessarily in equilibrium with daughters at time of shipment.

U-238D	AF-238D-A2	1,18 Bq	maximum activity, Not supplied as spectroscopy grade sources. Not necessarily in equilibrium with daughters at time of shipment.
--------	------------	---------	--

Back to overview - Alpha sources AF-Type

([home/products/isotope_products/isotrak_calibration_sources/reference_sources/planchet_disk_sources/alpha_sources_af_type/](#))

**USA / North and South America
/ Australia / China & Far East
Asia**

**Eckert & Ziegler
Analytics, Inc.
1380 Seaboard Industrial Blvd.
Atlanta, GA 30318 (USA)
Phone: +1 404 352 8677
Fax: +1 404 352 2837
E-Mail: isotrak.usa@ezag.com**

PRODUCT CONTACT

**Europe / Africa / Middle East
Asia / Japan**

**Eckert & Ziegler
Nuclitec GmbH
Gieselweg 1
38110 Braunschweig
Germany
Phone: +49 5307 932 0
Fax: +49 5307 932 272
E-Mail: infoisotrak@ezag.com**

PRODUCT CONTACT

HOW TO FIND US

([home/products/isotope_products/hc](#)

France

Eckert & Ziegler Isotope

Products France

12 avenue des Tropiques

Hightec Sud - bâtiment B

91955 Courtaboeuf cedex

France

Phone: +33 1 64 86 22 22

Fax: +33 1 69 86 10 65

E-Mail: rudym.operon@ezag.com

PRODUCT CONTACT

Copyright by Eckert & Ziegler. All rights reserved.

[IMPRINT \(HOME/IMPRINT/\)](#) [TERMS OF USE \(HOME/TERMS_OF_USE/\)](#) [DATA PROTECTION \(HOME/DATA_PROTECTION/\)](#) [SITEMAP \(HOME/SITEMAP/\)](#)
[CONTACT \(HOME/ABOUT_US/CONTACTFOLDER/CONTACT/\)](#)